



**MICROBS**

**CHECK FAST!**

**Your microbiology result in 15 minutes**



# IAN

Detects your microbiological contamination in just **15 minutes.**

To boost productivity in your lab, Microbs launched IAN, an **Ultra Rapid Microbiological Method** allowing an accurate enumeration within minutes, not days!

By combining benefits of **solid phase cytometry** and **Artificial Intelligence**, Microbs overtaken the limits of current techniques.

Adjusted on more than **40 million events**, IAN<sup>®</sup> detects bacteria, yeasts and molds without the need for incubation or enrichment.

Used at line testing or in your laboratory, this **fully automated technology** will prevent and reduce risks of contamination in your products.



Pharmaceutical



Cosmetics



Food & beverage

IAN



## Ultra-rapid microbiological method

- ✓ Bacteria, Yeasts and Molds detection
- ✓ Preprogrammed cap
- ✓ Time to results in 15 min.
- ✓ Fully automated
- ✓ No need for interpretation

### Increase your productivity:



Optimize  
your process



Decrease your  
product losses



Reduce your  
inventory




Increased  
efficiency

Check fast!





**Faster anywhere,  
for everyone.**  
**IAN, a breakthrough  
technology**



Touchscreen

LED status indicator strip

Carrying Handle

Analysis Buffer

Decontamination Buffer

USB port for data export

RFID detector for enhanced traceability



**Mobile Laboratory**

Autonomous, it fits where you need it most



**Simplicity**

No expertise required



**Speed**

No enrichment needed



**Accuracy**

Detection range from 1 to ≈ 5 000 microorganisms/volume analyzed



**Liquid sample**

Sampling from 100µL to 200mL



**Traceability**

Full Data integrity  
21 CFR 11 compliant

**IAN Dimensions:** Open drawer 24 x 46 x 52 cm. - 9,5" x 18,1" x 20,5"  
**Weight:** 16kgs - 35.3 pounds - **Power supply:** 110-230V.  
No IT required, can operate off network. Secured and proprietary system.

**MICROBS**  
6, rue Pierre Joseph Colin  
35000 RENNES - France

**+33-230-963-641**  
**www.microbs.fr**





## Pharmaceutical industry:

### Take control of all key steps in your production process

Microbs disrupts the field of bioburden and sterility testing.

To quickly measure microbial contamination in your product, Microbs launched IAN®, an innovative direct non growth counting method.

This unique association of Solid Phase Cytometry and Artificial Intelligence allows ultra rapid detection of aerobic and anaerobic microorganisms in your sample in only 15 min.

## DEPARTMENTS

### MANUFACTURING

#### Optimize your process

- Reduce lead time
- Improve production capacity



SAVE UP TO 14 DAYS ON QUARANTINE TIME

## APPLICATIONS

### Bioburden test

- Intermediate products
- Raw material
- WFI & Process water



### Sterility test

- Sterility of filterable final products



### Environmental control

- Swabbing surface / people
- Air monitoring



### Out of Specification Investigation

- Ultrafast investigation of sources of contamination



### QUALITY

#### Reduce Quality costs

- No defects
- No recalls



### LABORATORY

#### Faster batch release

- Rapidity
- Automation
- Traceability



### SUPPLY CHAIN

#### Reduce inventory

- Save inventory costs
- Reduce safety stock & required storage area



Increased product safety

Check fast!





# Microbiological testing is now fast and easy!

Just insert your caps with your sample and let IAN do the rest



TIME TO RESULT:  
**15 mn.**



## What happens inside:

- 1 **FILTRATION:** Sample aspiration through microfluidics for microorganisms retention.
- 2 **LABELLING:** Sequential implementation of unique formulation of dyes.
- 3 **DETECTION:** Data are collected all through the process by multiple sensors.
- 4 **ANALYSIS:** over 100 criteria analyzed by the AI, trained on 40 million events.
- 5 **RESULTS:** Quantification of only viable microorganisms. Certificate of analysis.

## IAN, a fully automated solution

### Accurate & validated method




Figure 1. Correlation between Microbs TVC Cap and Compendial method. The straight-line  $y=x$  is shown as a marker of linearity. The red points correspond to Pharmacopoeial strains.



### 165 Bacteria (including 65 genera)

Such as *Pseudomonas aeruginosa*, *Escherichia coli*, *Staphylococcus aureus*, *Bacillus subtilis*, *Clostridium sporogenes*



### 37 Yeasts & Molds (including 20 genera)

Such as *Aspergillus niger*, *Candida albicans*, *Saccharomyces cerevisiae*, *Candida parapsilosis* and *Zygosaccharomyces bailii*



Compliant with European (EP.5.1.6) & US pharmacopeia (USP 1223) validation:

Bioburden: Accuracy, Linearity, LOQ, Specificity and Robustness.  
Sterility test: LOD, Specificity and Robustness.



## Services

### Documentation & Execution

From feasibility testing to regulatory documentation, our dedicated professional services are the single point of contact to fit all your needs:

### Feasibility study

### IQ & OQ

### Method validation

(Qualitative & Quantitative)

### Product validation

(Qualitative & Quantitative)

MICROBS  
6, rue Pierre Joseph Colin  
35000 RENNES - FRANCE

+33-230-963-641  
www.microbs.fr

MICROBS



**MICROBS**

6, rue Pierre Joseph Colin — 35000 RENNES - FRANCE

☎ **+33-230-963-641** — [www.microbs.fr](http://www.microbs.fr)