



# SYNCT™

SYNCT™ Software | The Power Behind Your Lab

## SYNCT™ Software

SYNCT™ Software empowers you to manage data from any place in the lab. It is designed specifically for molecular labs running ARIES® and NxTAG® assays, and provides flexible, customizable workflows by connecting laboratory workstations. SYNCT offers:

### Connectivity


- Use a bidirectional interface to connect SYNCT-enabled Luminex Systems to your laboratory information system (LIS) using HL7 or CSV
- Integrate additional instruments to your LIS in just a few clicks

### Data Management

- View reports from anywhere in your lab
- Create and prepare orders away from the instrument
- Obtain data for multiple tests in one place with centralized reporting

### Flexibility

- Adapt workflow to fit your lab process
- Expand functionality to your software with SYNCT applications (apps)
- Expand your menu by running LDTs developed using the User Defined Protocol (UDP) app on SYNCT



## Recommended Hardware Specifications

### SYNCT Software (without xPONENT®)

- Microsoft Windows® 7 Professional, SP1 (32-bit or 64-bit, US English), or Microsoft Windows 8.1 Update (64-bit, US English)
- CPU: 2 gigahertz (GHz) or faster
- 4 GB or greater RAM
- 160 GB hard drive (or higher)
- Graphic device with WDDM drive at minimum 1280 x 1024 resolution

### SYNCT Software (with xPONENT)

- Microsoft Windows 7 Professional, SP1 (64-bit, US English)
- CPU: 2 gigahertz (GHz) or faster
- 8 GB or greater RAM
- 160 GB hard drive (or higher)
- Graphic device with WDDM drive at minimum 1280 x 1024 resolution

## Ordering Information

Description	Part Number
SYNCT™ (RUO)*	CN-SW34
SYNCT™ (IVD)**	CN-SW47

\*For Research Use Only. Not for use in diagnostic procedures.

\*\* For In Vitro Diagnostic Use. Products are region specific and may not be approved in some countries/regions.

**Luminex**  
complexity simplified.

**orders@luminexcorp.com or support@luminexcorp.com**

Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. Validation of the LIS compatibility must be performed by the end user. ARIES® Systems are class 1(I) laser products.

©2017 Luminex Corporation. All rights reserved. Luminex, ARIES, NxTAG, and xPONENT are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

[www.luminexcorp.com/synct](http://www.luminexcorp.com/synct)

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europe@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.802.31150

info@luminexcorp.com

JAPAN

+81.3.5545.7440

infojp@luminexcorp.com

SS19821