

3M Food Safety

3M™ Sample Collection, Preparation and Processing.

Consistency from the start.

The complete 3M sample handling system.



3M[™]Sample Collection, Preparation and Processing products.

Better testing begins with better samples – and better samples are obtained when you use the best collection tools. 3M™ Sample Collection, Preparation and Processing solutions provide quality products for every step of your process. Convenient, ready-to-use products. Consistent, reliable results. Every time.











Biocides



Environment



Strength



Consistency



Temperature

Featured 3M[™] Sample Collection, Preparation and Processing products.





- Stick-mounted sponge
- Provides the ability to easily snap off the sponge from the handle after swabbing
- Available in wet or dry format
- ▶ 1.5" x 3" Biocide-free sponge



3M™ Hydrated-Sticks

- ▶ 1.5" x 3" biocide free cellulose sponge
- Available with or w/out gloves
- Pre-hydrated with range of media
- ► Bag with write-on-label
- ▶ Decontaminated through gamma irradiation
- Dry-sponge available



3M[™] Sample Bags

- Available with or w/out filter
- Bag with write-on-label
- ► Tear-off opening remains sealed until ready for use









3M[™] Flip-Top bottles 3M[™] Mini Flip-Top vials

- Volume 90 ml/bottle, Mini Flip-Top 9 ml/vial
- Hinged lid stays open for easy sample additions
- Ready to use with no additional blending required
- Headspace and bottle opening make sample addition easy
- Leak resistant vial helps prevent leakage during shipping
- Available with various diluents

3M[™] Swab-Samplers

- Swab reaches into corners and hard-to-reach areas for increased accuracy
- Convenient design helps improve sampling consistency from tech-to-tech and plant-to-plant
- Bottles have ample headspace for sample addition and vortexing
- Self-contained swab and vial for increased productivity and efficiency
- Swab in pre-filled bottle with various diluents

3M™ Quick Swab

- Ready to use environmental swab quickly delivers 1 ml of sample onto a 3M™ Petrifilm™ Plate
- Self-contained sampling swab is easy to use and can reduce need for dilution mixing and pipetting for quick sampling
- Contains Letheen neutralising buffer to facilitate recovery of bacteria

3M™ Enviro Swab

- Ready to use swab reaches into corners and hard-to-reach areas for increased accuracy
- ► Temper-evident band indicates if the swab has been used
- Convenient design helps improve sampling consistency from tech-to-tech and plant-to-plant
- Resealabe tube provides tight closure for sample transportation



Transport.

Wide variety of bags ensures collected samples are uncompromised during transport to testing lab.

Sample bags

- ► Durable, quality construction
- ▶ Wire ties for strong seals and reduced leakage
- ► Easy-to-open bag
- Available with or without sponges

Prepare.

Prepared enrichments and diluents for standardised results – test after test.

Prepared enrichment media

- Cost-effective alternative to making and sterilising media
- ► Lowers operational labour costs
- ► Reduces errors in mixing and preparing dilutions
- Wide range of broths and buffers





Analyse.

The comprehensive line of 3M Food Safety products works together to deliver consistent, reliable results.

Testing solutions and lab supplies

- ➤ 3M[™] Petrifilm[™] Plates Sample-ready microbiology indicator count plates
- ► 3M[™] Molecular Detection System
 Pathogen testing with molecular accuracy
- **► 3M[™] Allergen Testing**



Sponges	3M™ Hydra-sponge	
	HS10NB2G	3M™ Hydra-Sponge Neut Buffer w/gloves, 5 bags of 20
	HS10BPW	3M™ Hydra-Sponge w/buffered Peptone Water Broth, 10ml
	HS10NB	3M™ Hydra-Sponge 10ml Neut Buffer 5BG OF 20, 100/CA
	HS10BPW2G	3M™ Hydra-Sponge/Buffpeptone Water Broth Gloves, 10ml
	HS10LET2G	3M™ Hydra-Sponge w/letheen Broth w/gloves, 10ml,
	3M™ Sponge stick	
	SSL10NB	3M™ Sponge-Stick w/neutralizing Buffer, 5 bags of 20
	SSL10BPW	3M™ Sponge-Stick w/buffered Peptone Water Broth, 10ml
	SSL100	3M [™] Sponge-Stick , 100/CA
	BP133ES	3M [™] Dry-Sponge, 18oz, 100/CA
	SSL10LET	3M [™] Sponge-Stick w/letheen Broth, 10ml, 100/CA
	3M™ Dry sponge	
	BP133ES	3M [™] Dry-Sponge, 18oz, 100/CA
Swabs	3M™ Enviroswab	
	ENVSWB100	3M™ ENVSWB100 Envrio Swab 100/CS
	3M™ Quick swab	
	6433	3M™ Quick Swab, 250 EA/Case
	Other swab sampler	
	RS9601LET	3M™ Swab-Sampler with Letheen Broth, 1ml, 100/CA
	RS96010BPW	3M [™] Swab-Sampler w/buffered Peptone Water Broth, 10ml
	RS96010LET	3M [™] Swab-Sampler with Letheen Broth, 10ml, 100/CA
	RS9604LET	3M [™] Swab-Sampler with Letheen Broth, 4ml, 100/CA
	RS9604NB	3M™ Swab-Sampler w/neutralizing Buffer, 4ml, 100/CA
	RS96010NB	3M™ Swab-Sampler w/neutralizing Buffer, 10ml, 100/CA
Dilutions	3M™ Flip top	
	BPPFV9BPW	3M™ Mini Flip Top Vial w/buffered Peptone Water Broth
	BPPFV9MRD	3M™ Mini Flip Top Vial w/max Recovery Diluent
	BPPFV9BFD	3M™ Mini Flip Top Vial w/Butterfield Buffer, 100/CA
	FTPW9060	3M [™] Flip Top Dilution Bottles w/Peptone Water Diluent, 60/CA
	FTBFD9060	3M™ Flip Top Dilution Bottles w/Butterfield Buffer, 60/CA
	FTBPW9060	3M™ Custom Flip Top Vial w/buffered Peptone Water Broth
	Enrichment pouch	
	QEBPWI225F	3M™ ENR PCH 225BPW ISO FTR 32CS
Sample bags	BP113S	3M™ Sample Bags, 22oz, 5.5x9", plain, 1000/CS
	6469	3M [™] Filter Bag, 200/box
	6468	3M™ Sample Bag, 500/box
	BP701	3M™ Sample Bags 55oz, 7x12", 3mil, 1000/CS

Global experience. Local support.

3M Food Safety products are at work in more than 60 countries around the globe, to help keep businesses like yours moving forward. Our passion is to work with you to discover new food safety solutions that protect your brand, mitigate your risk and improve your operations. And, every day, our people are there to support you locally, by providing technical help and exceptional service.

For more information, contact your local 3M Food Safety representative or visit **3M.com/foodsafety**





3M UK PLC Charnwood Campus 10 Bakewell Road Loughborough LE11 5RB

emeafoodsafety@mmm.com www.3M.com/foodsafety